
LiNA OperåScope™

Biopsy Forceps


Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

1


OP-FCP-6


2


3a


3b


3c


4a b


4c

4d

4e f

5

LiNA  LiNA OperâScope™ Biopsy Forceps

REF OP-FCP-6

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Product description:

The LiNA OperâScope™ Biopsy Forceps (Forceps) are delivered as a sterile, single-use device. The Forceps are designed to be used specifically with the LiNA OperâScope™ endoscope. The device consists of a handle, a shaft, and a jaw system (forceps) at the distal end. The Forceps are compatible with endoscopes with a straight working channel of min. Ø1.87 mm diameter and a max length of 310 mm.

Intended use:

LiNA OperâScope™ Biopsy Forceps are intended for use in hysteroscopy and female cystoscopy procedures to obtain endometrial or urinary bladder tissue samples for histologic examination.

Indications for use:

- Benign and malignant endometrial diseases, such as endometrial polyps, myoma or endometrial cancer.
- Urothelial abnormalities such as bladder cancer/carcinoma in situ and interstitial cystitis.

Patient population:

Adolescent and adult females suited for hysteroscopy and cystoscopy.

Contraindications:

This product is not intended for use when hysteroscopy or cystoscopy is contraindicated.

Contraindications for hysteroscopy include:

- Acute Pelvic Inflammatory Disease (PID)

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- Inability to distend the uterus,
- Cervical stenosis,
- Cervical/vaginal infection,
- Uterine bleeding or menses,
- Known pregnancy,
- Invasive carcinoma of the cervix,
- Recent uterine perforation,
- Medical contraindication or intolerance to anesthesia.

Endometrial Biopsy:

- Clotting disorder or coagulopathy

Cystoscopy:

- Acute/known Urinary Tract Infections,
- Severe Coagulopathy.
- Medical contraindication or intolerance to anesthesia.

Potential complications of hysteroscopy:

- Potential complications of Continuous Flow Hysteroscopy and Cystoscopy include: uterine or bladder perforation/false passage (resulting in possible injury to bowel, bladder, major blood vessels and ureter), cervical laceration, hemorrhage, vasovagal episodes, pain, infection, uterine adhesions, air embolism and fluid-related complications (hyponatremia, hypothermia, pulmonary edema, cerebral edema, heart failure), urethral trauma, irritable bladder symptom. Endometrial biopsies can lead to cramping and vaginal bleeding/spotting.

Warnings:

- The LiNA OperâScope™ Biopsy Forceps are provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- For single use only. Do not reuse, reprocess or re-sterilize the LiNA OperâScope™ Biopsy Forceps. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device result in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA OperâScope™ Biopsy Forceps.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA OperâScope™ Biopsy Forceps should

LiNA  LiNA OperâScope™ Biopsy Forceps

REF OP-FCP-6

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

not be used for any other purpose than their intended function.

- To mitigate the risk of perforation, only advance or manipulate the LiNA OperâScope™ Biopsy Forceps while viewing a LiNA OperâScope™ live camera image, allowing observation of the uterine cavity or urinary bladder.
- Do not advance the device should resistance be experienced while in the patient.
- Consult the Instructions for Use of the endoscope before performing procedures with the LiNA OperâScope™ Biopsy Forceps.

Precautions:

- Use of the LiNA OperâScope™ Biopsy Forceps requires adequate training and experience in performing hysteroscopy or cystoscopy.
- Always have a backup device readily available for immediate use. If any malfunction should occur during use, stop the procedure immediately, and slowly withdraw the LiNA OperâScope™ Biopsy Forceps and replace with a new device.
- After use, the device LiNA OperâScope™ Biopsy Forceps is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances (**Figure 5**).
- The LiNA OperâScope™ Biopsy Forceps are intended only as an adjunct in assessing patient condition. It must be used in conjunction with clinical signs and symptoms.

Instruction for use:

- Inspect the label to ensure that the correct device is chosen (**Figure 1**).
- Using sterile technique, remove the device from the sterilized blister (**Figure 2**).
- Remove the protection cap (**Figure 3a**).
- Inspect the Forceps shaft and jaws for any obvious damage.
- Check the open/close function of the Forceps by actuating the handle. Opening and closing the Forceps is done by pushing or pulling the Handle Knob away from/ towards the Finger Eyelet (**Figure 3b & 3c**).

6. Close the Forceps (**Figure 4a**) before inserting the device into the LiNA OperâScope™ working channel (**Figure 4b**).

7. Slowly insert the Forceps into the LiNA OperâScope™ working channel until the jaws are visualized on the LiNA OperâScope™ live camera image. Some resistance will be felt when the jaw system passes through the LiNA OperâScope™ pre-curved tip. Take care not to bend the Forceps shaft.

8. Maneuver the Forceps towards the targeted site. Open the Forceps by activating the handle (**Figure 4c**). Advance the opened Forceps over to the targeted site, rotate the Forceps, if needed, (**Figure 4d**) and close the Forceps. Use only enough force to bite the tissue. Over-exertion could cause the Forceps to fail.

9. Continue to apply gentle force on the handle (**Figure 4e**) as the Forceps are extracted from the working channel (**Figure 4f**). If the tissue sample is too large to be extracted through the working channel, remove the Forceps and the LiNA OperâScope™ simultaneously.

10. Remove and prepare the tissue sample according to standard technique for histologic evaluation.

11. Steps 4 – 10 may be repeated if multiple samples are needed within a single procedure.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/ or patient is established.

Country of origin:

Poland.

Original IFU
2025-10 FV0254F

Original IFU
2025-10 FV0254F

 LiNA Medical AG - Platz 3 - 6039 Root D4 - Switzerland - Tel: +41 78 617 56 56 - e-mail: info@lina-medical.com - www.linamed.com
LiNA is a registered trademark of LiNA Medical in EU, US and CN

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