LiNA OperåScope™

Rat Tooth Alligator Grasper

1

OP-GRA-6

2

LINA

3 a

3 b

3 c

4 a b

view

4 c

view

4 d

view

4 e f

view

5

Recycling symbol

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Product description:

The LiNA OperaScope™ Rat Tooth Alligator Grasper (Grasper) is delivered as a sterile, single-use device. The LiNA OperaScope™ Rat Tooth Alligator Grasper is designed to be used specifically with the LiNA OperaScope™ endoscope. The device consists of a handle, a shaft, and a jaw system (rat tooth alligator grasper) at the distal end. The Grasper is compatible with endoscopes with a straight working channel of min. Ø1.87 mm diameter and a max length of 310 mm.

Intended use:

The LiNA OperaScope™ Rat Tooth Alligator Grasper is intended to be used in hysteroscopy procedures for retrieval of foreign bodies and endometrial polyps from the uterus and in female cystoscopy procedures for retrieval of foreign bodies and bladder stones from the urinary bladder.

Indications for use:

- Any need to retrieve endometrial polyps or foreign bodies from the uterus during hysteroscopy procedures.
- Any need to retrieve bladder stones or foreign bodies from the urinary bladder during female cystoscopy procedures.

Patient population:

Adolescent and adult females suited for hysteroscopy and cystoscopy.

Contraindications:

Contraindications for hysteroscopy include:

- Acute Pelvic Inflammatory Disease (PID)

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- Inability to distend the uterus,
- Cervical stenosis,
- Cervical/vaginal infection,
- Uterine bleeding or menses,
- Known pregnancy,
- Invasive carcinoma of the cervix,
- Recent uterine perforation,
- Medical contraindication or intolerance to anesthesia.

Cystoscopy:

- Acute/known Urinary Tract Infections,
- Severe coagulopathy,
- Medical contraindication or intolerance to anesthesia.

Potential complications:

Potential complications of Continuous Flow Hysteroscopy and Cystoscopy include: uterine or bladder perforation/false passage (resulting in possible injury to bowel, bladder, major blood vessels and ureter), cervical laceration, hemorrhage, vasovagal episodes, pain, infection, uterine adhesions, air embolism and fluid-related complications (hyponatremia, hypothermia, pulmonary edema, cerebral edema, heart failure), urethral trauma, irritable bladder symptom.

Warnings:

- The LiNA OperaScope™ Rat Tooth Alligator Grasper is provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product's sterile barrier or its packaging is damaged.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- For single use only. Do not reuse, reprocess or re-sterilize the LiNA OperaScope™ Rat Tooth Alligator Grasper. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA OperaScope™ Rat Tooth Alligator Grasper.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA OperaScope™ Rat Tooth Alligator

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

Grasper should not be used for any other purpose than intended function.

- To mitigate the risk of perforation, only advance or manipulate the LiNA OperaScope™ Rat Tooth Alligator Grasper while viewing a LiNA OperaScope™ live camera image, allowing observation of the uterine cavity or urinary bladder.
- Do not advance the device should resistance be experienced while in the patient.
- Consult the Instructions for Use of the endoscope before performing procedures with the LiNA OperaScope™ Rat Tooth Alligator Grasper.

Precautions:

- Use of the LiNA OperaScope™ Rat Tooth Alligator Grasper requires adequate training and experience in performing hysteroscopy or cystoscopy.
- Always have a backup device readily available for immediate use. If any malfunction should occur during use, stop the procedure immediately, and slowly withdraw the LiNA OperaScope™ Rat Tooth Alligator Grasper and replace with a new device.
- After use, the device LiNA OperaScope™ Rat Tooth Alligator Graspers is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/ or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances (**Figure 5**).
- The LiNA OperaScope™ Rat Tooth Alligator Grasper is intended only as an adjunct in assessing patient condition. It must be used in conjunction with clinical signs and symptoms.

Instruction for use:

1. Inspect the label to ensure that the correct device is chosen (**Figure 1**).
2. Using sterile technique, remove the device from the sterilized blister (**Figure 2**).
3. Remove the protection cap (**Figure 3a**).
4. Inspect the Grasper shaft and jaws for any obvious damage.

5. Check the open/close function of the Grasper by actuating the handle. Opening and closing the Grasper is done by pushing or pulling the Handle Knob away from/towards the Finger Eyelet (**Figure 3b & 3c**).

6. Close the Grasper (**Figure 4a**) before inserting the device into the LiNA OperâScope™ working channel (**Figure 4b**).
7. Slowly insert the Grasper into the LiNA OperâScope™ working channel until the jaws are visualized on the LiNA OperâScope™ live camera image. Some resistance will be felt when the jaw system passes through the LiNA OperâScope™ pre-curved tip. Take care not to bend the Grasper shaft.
8. Maneuver the Grasper towards the targeted site. Open the Grasper by activating the handle (**Figure 4c**). Advance the opened Grasper over to the targeted site, rotate the Grasper, if needed, (**Figure 4d**) and close the Grasper. Use only enough force to grasp the polyp or object. Overexertion could cause the Grasper to fail.
9. Continue to apply gentle force on the handle (**Figure 4e**) as the Grasper is extracted from the working channel (**Figure 4f**). If the tissue or object is too large to be extracted through the operative channel, remove the Grasper and the LiNA OperâScope™ simultaneously.
10. Remove and prepare the tissue or object according to standard technique for histologic evaluation and/or foreign body retrieval.
11. Steps 4–10 may be repeated if multiple retrievals are needed within a single procedure.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

Country of origin: Poland.

er

LiNA OperâScope™ Rat Tooth Alligator Grasper

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