
LiNA OperåScope™

Lasso

1

OP-LAS10-3
OP-LAS16-3

2

3 a

3 b

3 c

4 a b

4 c

4 d

4 e f

5

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Product description:

The LiNA OperâScope™ Lasso (Lasso) is delivered as a sterile, single-use device. The Lasso is designed to be used specifically with the LiNA OperâScope™ hysteroscope. The device consists of a handle, a shaft and the basket at the distal end. The Lasso is available in two different polyp basket sizes, the OP-LAS10-3 with 10 mm diameter and OP-LAS16-3 with 16 mm diameter. The Lasso is compatible with hysteroscopes with a straight working channel of min. Ø1.87 mm diameter and a max length of 310 mm.

Intended use:

The LiNA OperāScope™ Lasso (Lasso) is intended for use in hysteroscopy procedures for retrieval of soft tissue such as polyps from the uterine cavity.

Indication for use:

Use in hysteroscopy procedures for retrieval of soft tissue such as polyps from the uterine cavity.

Patient population:

Adolescent and adult females suited for hysteroscopy.

Contraindications:

This product is not intended for use when hysteroscopy is contraindicated. Contraindications for hysteroscopy include:

- Acute Pelvic Inflammatory Disease (PID)

Hysteroscopy may be contraindicated by the following conditions, depending on their severity or extent:

- Inability to distend the uterus,
- Cervical stenosis,
- Cervical/vaginal infection,
- Uterine bleeding or menses,
- Known pregnancy,
- Invasive carcinoma of the cervix,
- Recent uterine perforation,
- Medical contraindication or intolerance to anesthesia.

Potential complications:

Potential complications of Continuous Flow Hysteroscopy include: uterine perforation/false passage (resulting in possible injury to bowel, bladder, major blood vessels and ureter), cervical laceration, hemorrhage, vasovagal episodes, pain, infection, uterine adhesions, air embolism and fluid-related complications (hyponatremia, hypothermia, pulmonary edema, cerebral edema, heart failure).

Warnings:

- The LiNA OperâScope™ Lasso is provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- For single use only. Do not reuse, reprocess or re-sterilize the LiNA OperâScope™ Lasso. Any reprocessing may impede the functions of this device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device result in risk of device malfunction and/or erroneous pathology specimen collection due to residual tissue in the LiNA OperâScope™ Lasso.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA OperâScope™ Lasso should not be used for any other purpose than the intended function.
- To mitigate the risk of perforation, only advance or manipulate the LiNA OperâScope™ Lasso while viewing a LiNA OperâScope™ live camera image, allowing observation of the uterine cavity.
- Do not advance the device should resistance be experienced while in the patient.

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

- Consult the Instructions for Use of the hysteroscope before performing procedures with the LiNA OperåScope™ Lasso.

Precautions:

- Use of the LiNA OperāScope™ Lasso requires adequate training and experience in performing hysteroscopy.
- Always have a backup device readily available for immediate use. If any malfunction should occur during use, stop on the procedure immediately, and slowly withdraw the LiNA OperāScope™ Lasso and replace with a new device.
- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.
- The LiNA OperāScope™ Lasso is intended only as an adjunct in assessing patient condition. It must be used in conjunction with clinical signs and symptoms.

the tip of the shaft passes through LiNA OperaScope's pre-curved tip. Take care not to bend the basket shaft.

8. Maneuver the basket towards the targeted site. Open the basket by activating the handle (Figure 4c). Advance the opened basket over to the targeted tissue, rotate the basket, if needed, (Figure 4d) and close the basket. Use only enough force to grasp the object. Over-exertion could cause the basket to fail.
9. Continue to apply gentle force on the handle (Figure 4e) as the basket is extracted from the working channel (Figure 4f). If the tissue is too large to be extracted through the operative channel, remove the Lasso and the LiNA Operascope™ together.
10. Remove and prepare the tissue according to standard technique for histologic evaluation.
11. Steps 4-10 may be repeated if multiple retrievals are needed within a single procedure.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

Country of origin: Poland.

Instructions for use:

1. Inspect the label to ensure that the correct device is chosen (Figure 1).
2. Using sterile technique, remove the device from the sterilized blister (Figure 2).
3. Remove the protection cap (Figure 3a).
4. Inspect the basket shaft and basket for any obvious damage.
5. Check the open/close function of the basket by actuating the handle (Figure 3b & 3c). Opening and closing the device is done by pushing or pulling the Handle Knob away from/ towards the Finger Eyelet.
6. Close the basket (Figure 4a) before inserting the device into the LiNA OperâScope™ working channel (Figure 4b).
7. Slowly insert the basket into the LiNA OperâScope™ working channel until the basket is visualized on the LiNA OperâScope™ live camera image. Some resistance will be felt when

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