




LiNA SeaStar™

(Sterile)

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

LiNA SeaStar™

REF 4000-S-SL, 4000-S-VL, 4000-S-VS,
4000-S-SS, 4000-S-VV, 4000-S-Q

en

Product description:

The LiNA SeaStar™ is a retractor ring with grooves every 2.5 cm on the convex site. The hinges are controlled by finger screws allowing the two ring halves to remain at any desired angle. The device is provided in a broad range of shapes and sizes which can easily be adjusted to fit various anatomical areas, thereby facilitating different surgical applications. The LiNA SeaStar™, together with LiNA Skin Hook™ is a complete system for retracting delicate soft tissue. The flexible design of the system allows easy adjustment and repositioning of the LiNA SeaStar™ during the surgical procedures and provides optimal exposure of the surgical site.

Intended use:

The LiNA SeaStar™ is intended to retract soft tissue from the operative surface.

Indication for use:

The LiNA SeaStar™ is intended to retract soft tissue from the operative surface.

Patient population:

All groups dedicated to surgical procedures which require optimal exposure.

Warnings:

- The LiNA SeaStar™ is supplied sterile using an ethylene oxide sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- For single use only. Do not reuse, reprocess or re-sterilize the LiNA SeaStar™. Any reprocessing may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- No modification of this equipment is allowed.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA SeaStar™ should not be used for any other purpose than intended function.

Precautions:

- Do not use excessive force during the screw tightening.

- Use of the LiNA SeaStar™ requires adequate training and experience in performing surgical procedures. Intended user profile is a trained surgeon or physician.
- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.
- Always have a backup device readily available for immediate use.

Instruction for Use:

Place the LiNA SeaStar™ around the initial incision and insert the hooks in directly opposing quadrants with equal tension. As dissection progresses, move the position of the hook into the deeper layers of the fascia. You may then add new hooks or reposition those already in place to achieve full balanced retraction from any angle. To ease closing maintain constant opposing tension.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

Caution:

Federal Law (USA) restricts this device to sale by or on the order of a trained physician.

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LiNA is a registered trademark of LiNA Medical in EU, US and CN