LiNA Gold Loop™

Innovation in Gynecology

LiNA Medical AG - Platz 3 - 6039 Root D4 - Switzerland - Tel: +41 78 617 56 56 - e-mail: info@lina-medical.com - www.linamed.com

LiNA Medical Polska Sp. z o.o. - Rolna 8A - Sady - 62-080 Tamowo Podgórze - Poland - Tel: +48 61 222 21 21 - e-mail: info.poland@lina-medical.com - www.lina-medical.pl

2025-09 FV0023R











LiNA is a registered trademark of LiNA Medical in EU, US and CN

2025-09 FV0023R

R Only

Refer to the Symbol Definition Chart FV0223, the integral part of the Instructions for Use, for the meanings of the symbols.

LiNA Gold Loop™

EL-160-4, EL-160-8, EL-200-4, EL-200-8, EL-240-4, EL-240-8

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

Product description:

The LiNA Gold Loop™ is a 5 mm monopolar electrosurgical device. It consists of wire loop inside an introducer sheath which can be inserted into the abdominal cavity through a laparoscopic access port.

The loop opens when the handle is pushed forward. Rhomboid shaped loop is available in different sizes depending on LiNA Gold Loop™ variant.

Intended use:

The LiNA Gold Loop™ is a monopolar device intended for Laparoscopic Supracervical Hysterectomy (LSH).

Indication for use:

The LiNA Gold Loop™ is a 5 mm monopolar electrosurgical device intended for cutting and removal of soft tissue in gynecological procedures involving endoscopic hysterectomy. It is used with standard electrosurgical generators.

Patient population:

Intended patient population are women undergoing Laparoscopic Supracervical Hysterectomy (LSH).

Contraindications:

The LiNA Gold Loop™ should not be used in the following situations:

- If the uterine isthmus cannot be visualized circumferentially before activation.
- If the uterus is adhered to any adjacent organs.
- If the patient has large posterior uterine or cervical wall myomas.
- Removal of uterus which are larger than a gestational week of 20.
- Removal of uterus which is attached to the pelvis with large laterally fibroids.
- The LiNA Gold Loop™ should not be used in patients who have undergone caesarean delivery due to the risk of bladder injury in the event that the bladder cannot be completely dissected off of the lower anterior uterine segment.

Potential complications:

- Electrosurgical complications (e.g., insulation failure, electrocution);
- Conversion to laparotomy;
- Hemorrhage, bleeding complications;
- Implantable electrical device malfunction;
- Small or large bowel injury, including perforation;
- Urinary tract injury;
- Unintended burns (might lead to hemorrhage, small or large bowel injury including perforation, urinary tract injury, conversion to laparotomy);
- Pain;
- Infection.

Warnings:

- The LiNA Gold Loop™ is provided sterile using ethylene oxide sterilization. Carefully inspect the packaging for any

damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.

- For single use only. Do not reuse, reprocess or re-sterilize the LiNA Gold Loop™. Any reprocessing may impede the functions of the device. Reusing single use devices may also increase the risk of cross contamination. Attempts to clean the device results in risk of device malfunction.
- Carefully inspect the device before use. If any damage or defect is identified, dispose of the device.
- No modification of this equipment is allowed. Do not attempt to modify the original rhomboid design of the LiNA Gold Loop™.
- Do not use if the loop is bent, twisted, malformed, or does not unfold completely to the original rhomboid design when the handle is completely advanced.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.
- Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol), as explosion may occur.
- Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire.
- Do not allow the cable of the LiNA Gold Loop™ to be in contact with skin of the patient or surgeon during electrosurgical activation.
- The LiNA Gold Loop™ shall be used with monopolar electrosurgical units, with a maximum operating voltage not exceeding 2.5kV peak, and may only be used with a IEC 60601-1, IEC 60601-2-2 and IEC 60601-1-2 approved generator.
- Use the lowest possible settings on the associated electrosurgical unit capable of achieving the desired surgical effect, in order not to discharge unnecessary high energy in the patient.
- The LiNA Gold Loop™ must be large enough to pass around the corpus of the uterus. Then apply the loop around the cervix.
- The LiNA Gold Loop™ must not be pre-tested up against bodily tissue prior to any cutting process.
- Always verify that the active (i.e. un-insulated) portion of the electrode is not in contact with any organ other than the uterus (e.g. bowel or vessel) before activating the generator.
- Do not use excessive force while using LiNA Gold Loop™. This could damage the product or result in entrapment and mechanical sectioning of the uterus with the risk of excessive bleeding.
- Do not use the LiNA Gold Loop™ for slicing the uterus as alternative to morcellation.
- Activate the active electrode only when it is in contact with target tissues and limit the time length of high-voltage peaks in order to reduce the risk for capacitive coupling.
- Avoid direct contact between the active electrode and other metal surgical equipment during surgery as this can result in thermal damage to tissue. Always keep the active electrode in view during the procedure.

LiNA Gold Loop™

EL-160-4, EL-160-8, EL-200-4, EL-200-8, EL-240-4, EL-240-8

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

- Remove the uterus manipulator before activating the LiNA Gold Loop™. Only uterus manipulator made of non-conducting material (e.g. plastic) should be used with the LiNA Gold Loop™.
- Remove intra uterine device (IUD) from the uterus prior to activating the LiNA Gold Loop™
- When not using LiNA Gold Loop™, place it in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.
- Due to concerns about the carcinogenic and infectious potential of electrosurgical byproducts (such as tissue smoke plume and aerosols), protective eyewear, filtration masks, and effective smoke evacuation equipment should be used.
- Aspirate fluid from the area before activating the instrument. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to an active electrode may carry electrical current or heat away from target tissues, which may cause unintended burns to the patient.
- Prior to increasing the intensity, check the adherence of the neutral electrode and its connections. Apparent low output or failure of the device to function correctly at the normal operating settings may indicate faulty application of the neutral electrode or poor contact in its connections.
- Pay attention when the patient has any implantable electronic device (IED). The perioperative team shall consult with the team managing the implanted device before surgery and as needed. In case team managing the device is not available, consult with the device manufacturer or interdisciplinary team in your hospital. When LiNA Gold Loop™ is activated it shall be kept more than ½ inch (1 cm) from the IED.
- If the patient has a metal implant, this must not be on the path between the neutral electrode and the LiNA Gold Loop™.
- Failure to follow all instructions or any warnings or precautions may result in significant injury to the patient, physician or attendants and may have an adverse effect on the outcome of procedures performed. LiNA Gold Loop™ should not be used for any other purpose than intended function.
- Connect the LiNA Gold Loop™ to the electrosurgical unit only when the energy is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.

Precautions:

- Always have a backup device readily available for immediate use. If any malfunction should occur during use, stop the procedure immediately, and slowly withdraw LiNA Gold Loop™ and replace with a new device.
- Use of the LiNA Gold Loop™ requires adequate training and experience in performing laparoscopic procedures.
- After use, the device is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.

Generator specifications and settings:

The LiNA Gold Loop™ is compatible with all standard electrosurgical generators with monopolar outlet and BLEND modes.

The recommended setting for the generator is 100 - 130 Watt, Blend 1. Max. voltage 2500 Vp-p.

Instruction for use:

After the uterine ligaments and vessels have been correctly ligated, when the uterus can be safely amputated from the cervix, follow the instructions below:

- Apply the neutral electrode to the patient and generator. Connect the LiNA Gold Loop™ to the generator's monopolar outlet. Set the generator at Blend 1 and start at 100 Watt.
- Carefully insert the LiNA Gold Loop™ through a 5 mm trocar/entry.
- Unfold the electrode (loop) by forwarding the handle. The handle must be completely advanced and the loop 100% unfolded prior to placing the loop around the uterine body.
 - Use grasper(s) and the uterine manipulator to help place the electrode (loop) around the uterine body.
 - Before activating the generator always verify that the active (i.e. un-insulated) portion of the electrode is not in contact with any organ other than the uterus (e.g. bowel or vessel).
- Carefully tighten the loop around lower part of uterine body (above cervix).
 - Remove the uterine manipulator and slightly release the tension of the loop.
- Start sectioning the uterine body by pressing "cut" on the foot pedal. After 1-3 seconds slowly pull the handle, until the electrode (loop) has cut the uterine body and is retracted inside the tube. Stop pressing the foot pedal.
- Remove the LiNA Gold Loop™ through the trocar.
- Dispose of the used LiNA Gold Loop™ after usage.
 - After amputating the uterus from the cervix, ensure that the proximal portion of the endocervical canal which may contain endometrial cells is coagulated, and ensure that any bleeding or oozing on the cervical stump is coagulated.
 - Provide instructions for follow up to patients in the event of complications including intra-abdominal bleeding which in rare cases can lead to exsanguinations.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

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LiNA is a registered trademark of LiNA Medical in EU, US and CH