
LiNA Xcise™

Laparoscopic Morcellator

Innovation in Gynecology

LiNA Medical AG - Platz 3 - 6039 Root D4 - Switzerland - Tel: +41 78 617 56 56 - e-mail: info@lina-medical.com - www.linamed.com


LiNA Medical Polska Sp. z o.o. - Rolna 8A - Sady - 62-080 Tarnowo Podgórne - Poland - Tel: +48 61 222 21 21 - e-mail: info.poland@lina-medical.com - www.lina-medical.pl

2025-10 FV0291D

LiNA is a registered trademark of LiNA Medical in EU, US and CN

2025-10 FV0291D

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

Product description:

The LiNA Xcise™ is a cordless, single use, fully disposable morcellator intended for tissue morcellation during laparoscopic gynecologic procedures. The LiNA Xcise™ contains a rotating cutting tube with built-in trocar function that also serves to protect the sharp end of the cutting tube. A grasper or a tenaculum forceps must be used to pull the stripes of tissue out through the lumen of the cutting tube. The LiNA Xcise™ cutting function is controlled by the activation button on the hand piece.

Intended use:

LiNA Xcise™ laparoscopic morcellator is intended for gynecologic use by trained professionals in hospital and surgical environments.

Indications for use:

Indicated for cutting, coring and extracting tissue in laparoscopic gynecologic procedures such as hysterectomy and myomectomy.

Patient population:

Females requiring laparoscopic gynecologic procedures such as hysterectomy and myomectomy. Should not be used in children.

Contraindications:

- Contraindications for use on vascularised tissue (ovaries, fallopian tubes, myomas and other structures): tissue must be devascularised and dissected before morcellation.
- Laparoscopic power morcellators are contraindicated in gynecologic surgery in which the tissue to be morcellated is known or suspected to contain malignancy.
- Laparoscopic power morcellators are contraindicated for removal of uterine tissue containing suspected fibroids in patients who are:
 - post-menopausal or over 50 years of age, or
 - candidates for en bloc tissue removal through the vagina or via a mini-laparotomy incision.

Potential complications:

- Direct trauma with injuries to surrounding tissue such as the bladder, ureter, bowel and vessels, haemorrhage, port-site hernia.
- Dissemination of benign tissue fragments, causing iatrogenic endometriosis, adenomyosis, leiomyomatosis, parasitic leiomyoma and leiomyomatosis peritonealis disseminate.
- Dissemination of malignant tissue fragments.

Warnings:

Uterine tissue may contain unsuspected cancer. The use of laparoscopic power morcellators during fibroid surgery may spread cancer, and decrease the long-term survival of patients. This information should be shared with patients when considering surgery with the use of these devices.

- The LiNA Xcise™ is provided sterile using radiation sterilization. Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged.
- Do not use past expiration date.
- Do not use if the device is exposed to non-sterile surfaces before procedure.

Refer to the Symbol Definition Chart FV0223, the integral part of Instructions for Use, for the meanings of the symbols.

en

- The use of LiNA Xcise™ does not pose any significant risk of reciprocal or electromagnetic interference. LiNA Xcise™ RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. If interference is suspected, move the interfering equipment away or increase the distance between the devices.
- Do not use excessive axial force while changing the trocar position from safety guard to Cut 1 and/or Cut 2 and remember to hold the bracket prior to switching mode as not following these may result in detaching trocar from the rest of the device.
- Be aware that the cutting tip of the LiNA Xcise™ is not in contact with other instrument for example grasping forceps aiming to hold the tissue in place while morcellating the target tissue. It may cause knife dulling/chipping.
- Do not use excessive force while using LiNA Xcise™. This could damage the product.
- Please be aware that reducer cap is only for inserting 5 mm instruments, do not use obturator when the reducer cap is mounted and do not push the cutting tube through the reducer cap. It will cause damage of the reducer cap.
- After use, the device and the battery is to be considered contaminated medical waste and may pose a potential infection risk. Be aware of sharp edges. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.

- To activate the cutting blade and begin morcellating, press the activation button on the hand piece while pulling pieces of tissue through the cutting tube.
- Release the activation button as soon as the stripe of tissue is extracted from the LiNA Xcise™.
- For use with a 5 mm instrument: Mount the reducer cap onto the back of the morcellator with a pressing and turning motion.
- After surgery, remove the LiNA Xcise™ from the abdominal cavity. For proper disposal, turn the trocar into the "Safety Guard" position, remove the reducer cap and insert the obturator. Open the bottom housing, by inserting a blunt instrument (e.g. clamp or tweezers) into the housing catch and twisting the instrument. Remove the battery pack. The device and the battery is to be considered contaminated medical waste and may pose a potential infection risk. Dispose of the device and packaging in accordance with applicable laws, administrative, and/or local government policy and hospital procedures, including those regarding biohazards, microbiological hazards and infectious substances.

Reporting:

Any serious incident that has occurred in relation to the device should be reported to the LiNA Medical AG and the competent regulatory authority of the country in which the user and/or patient is established.

Country of origin: Poland.

IP class: IPX0.

Instruction for use:

The surgeon should read the Instructions for use carefully before using this device.

- Carefully inspect the packaging for any damage prior to use. Do not attempt to use the device if product sterile barrier or its packaging is damaged. Prior to removing the device from the tray take off the retainer.
- Prior to using the trocar function of the LiNA Xcise™, insert the obturator fully into the device.
- Be sure that the trocar is placed in the "Safety Guard" position. If not, place the trocar in the "Safety Guard" position by holding in the bracket and then simultaneously turn the trocar (see illustration in picture 2).
- The LiNA Xcise™ with obturator should be placed into the abdomen using standard technique for laparoscopic trocar placement. It is recommended to insert the LiNA Xcise™ with obturator through a 12-14mm incision under direct visualization.
- In order to remove tissue use a 10-12 mm forceps or similar instrument inserted through the lumen of the LiNA Xcise™ and into the abdomen. To prevent injury to the abdominal wall, the tissue to be morcellated should be completely exposed before attempting to extract it through the morcellator.
- It is recommended to use a second pair of grasping forceps to hold the tissue in place and reduce tissue movement during morcellation.
- Place the trocar in the required position by turning the trocar into cutting position: "Cut 1" for peeling function or "Cut 2" for coring function. Prior to changing the required position hold the bracket and simultaneously turn the trocar.
- Adjust the coreguard if needed.

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