

Cystoscopy without  
capital cost or  
complexity



## LiNA OperåScope®

Single-use Cystoscopy System

The complete, portable system  
for cystoscopy.



# LiNA OperåScope®

Single-use Cystoscopy System

The LiNA OperåScope® single-use cystoscopy system is designed for use in visualization of the female urethra and the bladder cavity during diagnostic and therapeutic cystoscopic procedures.

- **Convenient solution** without capital cost or complexity
- **Single-use** design is ready for use out of the box
- **360° steerable, deflected tip** allows for precise control and maneuverability
- **Slim 4.2 mm catheter** designed for improved patient tolerability
- **5.5 french operative channel and single-use instrumentation** for therapeutic use

## Ordering information:

Product Code	Description	UOM
OP-201-6	LiNA OperåScope® with HDMI cable & LCD (sterile, single-use)	Box of 6ea
OP-FCP-6	LiNA OperåScope® Biopsy Forceps (sterile, single-use)	Box of 6ea
OP-GRA-6	LiNA OperåScope® Rat Tooth Alligator Grasper (sterile, single-use)	Box of 6ea
OP-NEE4-6	LiNA OperåScope® Needle, 4 mm (sterile, single-use)	Box of 6ea
OP-NEE8-6	LiNA OperåScope® Needle, 8 mm (sterile, single-use)	Box of 6ea
OP-RM-1	LiNA OperåScope® Recording Module	Each
OP-TDK-6	LiNA OperåScope® Tubing & Drape Kit (sterile, single-use)	Box of 6ea

 LiNA Medical AG - Platz 3 - 6039 Root D4 - Switzerland  
Phone: +41 78 617 56 56, e-mail: [info@lina-medical.com](mailto:info@lina-medical.com), [www.linamed.com](http://www.linamed.com)

©2025 LiNA Medical AG. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions. LiNA is a registered trademark of LiNA Medical in EU, US and CN.

**Distributor:**  
LiNA Medical USA  
Norcross, GA 30093, USA  
Tel: +1 (855) 546-2633  
Fax: +1 (678) 380-0500  
[info@linamed.com](mailto:info@linamed.com)  
[www.linamed.com](http://www.linamed.com)

EN ISO 13485:2016  
Quality Management System



Moustafa S, Rosen E, Goodman L. Patient and provider satisfaction with saline ultrasound versus office hysteroscopy for uterine cavity evaluation prior to in vitro fertilization: a randomized controlled trial. J Assist Reprod Genet. 2021 Mar;38(3):627-634.