

# LiNA OperåScope® Needle

*Single-use rigid needle, ideal for in-office hysteroscopic injections*




Innovation in Gynecology

# LiNA OperåScope® Needle


Single-use rigid needle, ideal for in-office hysteroscopic injections

## Simplicity and Versatility

- Sterile, single-use
- Two sizes – to accommodate uterine wall thickness
- Protection cap
- Limit buckle to control needle exposure during introduction
- Luer lock for syringe connection



## Available LiNA OperåScope® Needle sizes



## Technical Specifications:

LiNA OperåScope® Needle	
Required channel inner diameter	5 Fr
Working length	313 mm
Needle width	23G
Needle length	4 or 8 mm
Tip shape	3-bevel tip

## Ordering information:

Ordering information:		LiNA OperåScope® Needle
Product code	Description	UOM
OP-201-6	LiNA OperåScope® with HDMI cable & LCD (sterile, single-use)	Box of 6 ea
OP-NEE4-6	LiNA OperåScope® Needle, 4 mm (sterile, single-use)	Box of 6 ea
OP-NEE8-6	LiNA OperåScope® Needle, 8 mm (sterile, single-use)	Box of 6 ea

©2025 LiNA Medical AG. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions. LiNA is a registered trademark of LiNA Medical in EU, US and CN.

EN ISO 13485:2016  
Quality Management System

LiNA Medical AG, Platz 3, 6039 Root D4 | Switzerland  
Phone: +41 78 617 56 56, info@lina-medical.com | www.linamed.com

EU REP LiNA Medical Polska Sp. z o.o., Rolna 8A, Sady, 62-080 Tarnowo Podgórne | Poland  
Phone: +48 61 222 21 21, info.poland@lina-medical.com | www.lina-medical.pl

