


LiNA Xcise[®]

Cordless Laparoscopic Morcellator



Cut the Cord



Innovation in Gynecology

LiNA Xcise®

Cordless Laparoscopic Morcellator



Simplicity. Instant-use fully disposable design.

No assembly required. No drive box or cables.

Lightweight cordless design incorporates motor and battery in one device. Open the package, insert the obturator, and you are ready to start morcellating.


Safety. Control your blade exposure.

LiNA Xcise® features 3 blade positions to choose from to aid in safety and control during tissue morcellation.



Speed. Motor Technology.

High torque ~1000 RPM motor to power through tissue.



360° adjustable coreguard

Ergonomic handle with precision trigger
Less hand fatigue and improved control with the ability to rapidly pulse the trigger.




Snap-on 5mm Reducer Cap




LiNA Xcise®

Cordless Laparoscopic Morcellator



Inner Bevel - LiNA Xcise®

Sharpness. 300 series stainless steel blade features a larger inner bevel to prevent dulling from inadvertent contact with the tenaculum.



Outer bevel of other commercially available morcellators

Indications for use:

Indicated for cutting, coring and extracting tissue in laparoscopic gynecologic procedures such as hysterectomy and myomectomy.

Important:

Please refer to the Instruction for Use for complete instructions, contraindications, warnings and precautions.

Ordering Information:

Product code	Description	UOM
MOR-1515-6	LiNA Xcise® – Laparoscopic Morcellator 15mm, battery powered morcellator system with obturator. Sterile, single use.	Box of 6 ea

Distributor:

LiNA Medical USA

Norcross, GA 30093, USA

Tel: +1 (855) 546-2633 | Fax: +1 (678) 380-0500 | info@linamed.com | www.linamed.com

©2024 LiNA Medical AG. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions.

LiNA is a registered trademark of LiNA Medical in EU, US and CN.



LiNA Medical AG, Platz 3, 6039 Root D4 | Switzerland
Phone: +41 78 617 56 56 , info@lina-medical.com | www.linamed.com

EN ISO 13485:2016
Quality Management System

