

LiNA Xcise®

Cordless Laparoscopic Morcellator



Innovation in Gynecology

LiNA Xcise®

Cordless Laparoscopic Morcellator



Simplicity. Instant-use fully disposable design.

No assembly required. No drive box or cables.

Lightweight cordless design incorporates motor and battery in one device. Open the package, insert the obturator, and you are ready to start morcellating.




Safety. Control your blade exposure.

LiNA Xcise® features 3 blade positions to choose from to aid in safety and control during tissue morcellation.



Speed. Motor Technology.

High torque ~1000 RPM motor to power through tissue.



360° adjustable coreguard

Ergonomic handle with precision trigger
Less hand fatigue and improved control with the ability to rapidly pulse the trigger.


Removable battery to facilitate easy disposal once the procedure is completed.

Snap-on
5mm Reducer Cap



LiNA Xcise®

Cordless Laparoscopic Morcellator



Sharpness. 300 series stainless steel blade features a larger inner bevel to prevent dulling from inadvertent contact with the tenaculum.

Inner Bevel - LiNA Xcise®

Outer bevel of other commercially available morcellators



Removable battery. To ease disposal.

Once the procedure has been completed, the LiNA Xcise® has a removable battery to facilitate easy disposal to address environmental guidelines and policies.

Indications for use:

Indicated for cutting, coring and extracting tissue in laparoscopic gynecologic procedures such as hysterectomy and myomectomy.

Important:

Please refer to the Instruction for Use for complete instructions, contraindications, warnings and precautions.

Ordering information:

Product code	Description	UOM
MOR2-1515-6	LiNA Xcise® - Laparoscopic Morcellator 15mm, battery powered morcellator system with obturator. Sterile, single use.	Box of 6 ea

Distributor:

LiNA Medical USA
 Norcross, GA 30093, USA
 Tel: +1 (855) 546-2633 | Fax: +1 (678) 380-0500 | info@linamed.com | www.linamed.com

©2025 LiNA Medical AG. All rights reserved. All products may not be available for sale in all geographies. Please contact your local LiNA Medical representative for information about products available in your area. This information is intended for healthcare professionals. If you are a patient seeking more information, please consult your healthcare provider. Before using any LiNA product, always refer to the IFU included in our packaging for complete indications for use, contraindications, warnings, precautions, and instructions.

LiNA is a registered trademark of LiNA Medical in EU, US and CN.